

Meet the CRE-E10, an easy, self-fitting hearing aid for those with mild to moderate hearing loss.

### **Wearer Benefits**

- For adults 18+ with perceived mild to moderate hearing loss.
- Earbuds are designed for comfort, ease, and everyday wear.
- Patient Customization using the Sony | Hearing control app:
  - Easily set up and self-adjust their hearing aids.
  - Control volume or sound quality for results tailored to their listening environment.
- Analyzes and optimizes sound based on the patients' surroundings to improve the listening experience in noisy environments.
- Bluetooth compatible (iOS only). Sony OTC hearing aids connect via Bluetooth to stream audio or listen to music.
- Uses a Qi-compatible rechargeable battery for wireless charging and up to 26 hours of continuous use. Fully recharged after 4 hours.
- Connects easily to both iOS and Android devices.

# SONY

#### **Kit Contents**


- 2 hearing aids (L&R)
- Charging Case
- 4 sizes of ear tips, devices are pre-fit with small ear tip
- Care Tools (cleaning brush, cloth)
- User Guide, Safety and Maintenance Information
- QR code for quick start and how to download the Sony | Hearing Control app
- Warranty Card

# **Return Policy**

- The wearer has 45 days to return undamaged products for credit.
  Both hearing aids and the charger must be included for a valid return.
- The return is valid for 135 days from the date of sale to the patient, inclusive of the 45 days for patient return.
- Sony OTC hearing aids must be returned to the seller from which they were purchased.

## Warranty

• Sony OTC hearing aids have a 1-year warranty that includes repair, but does not include loss and damages.



#### **CRE-E10 Specifications**

Colors	Black
User Self-fitting	Yes
Types of ear-tips that can be used	XS, S, M, L Open (Vent) type / Close type Optional: Open (Vent) type
Battery	Rechargeable
Mic # on each	2
Continuous usage time	Up to 26 hr after a full charge Up to 23 hr incl. 2 hr audio streaming
Indicator LED on Hearing Aid	N/A
Indicator LED on Charger	Two for hearing aid status One for charger status
IPX	IP68
OSPL 90 (Peak)	118 dB SPL
HFA-FOG	39 dB
Hearing Profiles	6
Streaming	iOS only (MFI)
Connectivity with Smartphone	Bluetooth
Smartphone App	iOS, Android

Apple App Store is a servicemark of Apple Inc. iOS is a trademark of Apple Inc., registered in the U.S. and other countries. Android and Google Play are trademarks of Google Inc.



